

## ComFort Heated Mitts

Heavyweight, machine washable polyester shell. Removable inserts contain a rice and herb blend.

Mitts can be used with a paraffin treatment to prolong the heat or as a stand-alone moist-heat treatment.

### CARE INSTRUCTIONS:

Microwave 1 1/2 - 2 min. or warm in hot towel cabi. Spot clean, or machine wash cover in cold water, gentle dry low heat.

**CAUTION: TEST TEMP BEFORE USE.**

### WARNING:

Bags are not toys. To avoid the danger of suffocation, keep plastic liners away from children, infants, and pets. Do not use in cribs, beds, carriages, or playpens.

WR MEDICAL ELECTRONICS CO.  
1700 Gervais Avenue  
Maplewood, MN 55109  
651-604-8400 | toll Free 800-321-6387

[www.therabath.com](http://www.therabath.com)



Rev. 09/15.